

STIMEL - 03 WHITE PAPER SERIES

---

White Paper 2

# EMG Biofeedback and Closed-Loop Neuromuscular Rehabilitation

*The Clinical Foundations of Intention-Driven Therapy*

---

## Abstract

---

Neuromuscular recovery after stroke and other neurological injuries depends on restoring voluntary motor control through repeated activation of damaged neural pathways. Many patients exhibit early motor intent that produces low-amplitude electromyographic (EMG) activity even when visible movement cannot yet occur. Closed-loop neuromuscular rehabilitation systems address this gap by detecting voluntary neuromuscular signals and synchronizing electrical stimulation with the patient's effort. By linking motor intent, muscle activation, and sensory feedback, these systems create a structured therapeutic cycle that supports neuroplastic reorganization. This white paper outlines the clinical and engineering foundations of EMG-based closed-loop therapy and describes how platforms such as Stimel-03 integrate signal detection, visual biofeedback, and intention-synchronized stimulation into a patient-driven rehabilitation model.

## Clinical Challenge in Early-Stage Motor Recovery

---

Voluntary motor control often persists even when overt movement is absent. After stroke, weakened descending neural drive may still recruit small motor units, producing EMG signals that reflect genuine motor intent. However, these signals are typically too weak to create joint movement, limiting the patient's ability to practice functional tasks during early therapy.

Because motor relearning depends on repeated activation of sensorimotor pathways, the inability to convert early neuromuscular activation into meaningful movement slows recovery. Technologies capable of detecting and reinforcing these subtle signals enable patients to participate actively in rehabilitation earlier, transforming intention into measurable therapeutic engagement.

## Principles of EMG Biofeedback

---

Surface electromyography provides real-time visibility into muscle activation. Even without visible movement, motor unit action potentials generate detectable electrical activity that can be processed and displayed to the patient.

EMG biofeedback improves motor learning by helping patients recognize successful muscle activation, encouraging graded modulation of voluntary effort, and increasing engagement by translating effort into immediate visual feedback.

Over repeated practice, biofeedback supports refinement of voluntary activation strategies and strengthens the neural circuits responsible for motor control.



**Figure 1. EMG Biofeedback Signal Detection** Surface EMG electrodes detect electrical activity generated by muscle fibers during voluntary contraction. The signal is processed and displayed as visual biofeedback, allowing patients and clinicians to observe neuromuscular activation in real time.

## Closed-Loop Neuromuscular Rehabilitation


Closed-loop rehabilitation systems extend EMG biofeedback by linking physiological monitoring directly to therapeutic intervention. Instead of simply displaying muscle activity, these systems use EMG signals to trigger electrical stimulation that assists the intended movement.

The rehabilitation cycle can be summarized as:

*Voluntary Intent → EMG Detection → Timely Stimulation → Movement Success → Sensory Feedback → Neural Reinforcement*

This temporal pairing between voluntary neural drive and assisted movement reinforces sensorimotor pathways according to well-established principles of activity-dependent neuroplasticity.

Synchronizing stimulation with voluntary activation strengthens associative learning mechanisms that support reorganization of motor pathways within the central nervous system.



**Figure 2. EMG-Triggered Closed-Loop Stimulation** Voluntary motor intent produces EMG activity in the target muscle. When the detected signal exceeds a predefined activation threshold, electrical stimulation is triggered to assist the intended movement. This closed-loop process links motor intent, muscle contraction, and sensory feedback, reinforcing neural pathways responsible for voluntary movement.

## Engineering Foundations of EMG-Triggered Therapy

Detecting early voluntary activation signals presents engineering challenges because EMG amplitudes during weak contractions may be extremely small and easily obscured by electrical noise, motion artifacts, and variations in skin impedance.

Advanced closed-loop rehabilitation systems therefore employ several signal-processing strategies:

- Signal filtering pipelines optimized for low-amplitude motor unit potentials
- Artifact rejection algorithms that suppress motion-related noise
- Adaptive activation thresholds that adjust to patient-specific signal patterns
- Low-latency signal processing to preserve synchrony between voluntary intent and stimulation delivery

These capabilities ensure that stimulation is triggered reliably and at the correct phase of the patient's voluntary effort, allowing weak neuromuscular activation to translate into meaningful therapeutic movement.

## Stimel-03: Intention-Driven Rehabilitation Platform

Stimel-03 integrates EMG detection, visual biofeedback, and electrical stimulation into a unified neuromuscular rehabilitation system designed for intention-driven therapy.

Surface electrodes detect voluntary neuromuscular activation attempts from the target muscle group. EMG signals are processed in real time and displayed visually so that patients can observe their effort during therapy.

When the detected EMG activity exceeds a predefined activation threshold, electrical stimulation is delivered to assist the intended movement. This synchronized interaction between voluntary activation and stimulation allows patients to convert weak neuromuscular signals into functional movement.

Clinicians can adjust stimulation parameters and activation thresholds according to the patient's rehabilitation stage. As voluntary motor control improves, stimulation support can gradually decrease while voluntary activation becomes the dominant driver of movement.

## BioRhythmIQ Signal Processing Architecture

---

The BioRhythmIQ signal processing architecture used in Stimel-03 is designed to detect extremely weak neuromuscular activation signals that often appear during early stages of recovery.


The system isolates low-amplitude motor unit action potentials through a structured signal-processing pipeline:

*Acquisition → Signal Filtering → Artifact Suppression → Adaptive Threshold Detection → Trigger Logic → Stimulation Output*

During acquisition, EMG activity from the target muscle is captured through surface electrodes. The signal then passes through filtering stages that isolate physiologically relevant frequency components while suppressing environmental electrical noise. Artifact-suppression algorithms remove disturbances caused by motion, electrode impedance shifts, and external electrical interference.

An adaptive detection threshold evaluates the processed signal to determine when genuine voluntary activation occurs. Once this threshold is exceeded, trigger logic initiates stimulation that assists the intended movement while preserving the temporal relationship between neural intent and motor output.

By enabling reliable detection of early voluntary activation attempts, the system allows therapeutic stimulation to reinforce motor intent even before visible movement becomes possible, allowing patients to begin active rehabilitation earlier in the recovery process.



**Figure 4. BioRhythmIQ Signal Processing Flow** Acquisition → Filtering → Artifact suppression → Adaptive threshold detection → Trigger logic → Electrical stimulation. The architecture converts weak neuromuscular signals into reliable therapeutic triggers while maintaining synchronization between patient effort and assisted movement.

## Unified Electrode Architecture

Stimel-03 uses the same electrode configuration for both EMG sensing and electrical stimulation. This unified architecture simplifies electrode placement and reduces therapy setup complexity.

Using a single electrode interface ensures that sensing and stimulation occur at the same anatomical location, improving signal consistency while reducing placement variability between sessions. The result is a streamlined clinical workflow with fewer technical steps for therapists.

## Clinical Applications and Use Cases

- Closed-loop EMG-triggered rehabilitation supports a range of neuromuscular recovery scenarios.
- Upper limb rehabilitation often targets wrist extension, finger opening, and grasp control. EMG-guided stimulation assists patients in performing functional reaching and grasping movements that are essential for daily activities.
- Lower limb rehabilitation commonly focuses on ankle dorsiflexion and knee stabilization during gait training. EMG-triggered stimulation can assist stepping movements and improve walking mechanics during therapy sessions.
- Closed-loop systems are particularly valuable during early recovery phases when EMG activity is present but visible movement is minimal.

## Clinical Impact

Synchronizing electrical stimulation with voluntary motor intent enhances patient engagement, increases repetition of meaningful motor activity, and reinforces sensorimotor pathways responsible for movement recovery.

By transforming weak voluntary activation into successful movement attempts, EMG-guided therapy allows patients to actively participate in generating therapeutic movement rather than receiving passive stimulation.

This approach supports progressive strengthening of voluntary motor pathways and contributes to functional recovery over time.

---

## Key Takeaways for Clinicians

- EMG-triggered functional electrical stimulation improves voluntary upper-limb motor recovery compared with passive stimulation approaches [3,4,8].
- Timing-synchronized stimulation enhances sensorimotor integration and promotes cortical reorganization following neurological injury [1,2].
- Detection of weak EMG signals allows patients with minimal visible movement to begin active rehabilitation earlier in the recovery process [4,6].
- Closed-loop neuromuscular training increases meaningful repetitions of motor activity, which is a key driver of neuroplastic adaptation and functional recovery [2,7].

---

## Conclusion

Closed-loop neuromuscular rehabilitation systems do more than assist movement. By pairing voluntary intent with assisted motion and sensory feedback, they create a continuous interaction between motor commands and sensory input. This interaction is fundamental to motor relearning because the nervous system depends on sensory information to refine and stabilize newly emerging motor patterns.

Understanding how biofeedback and proprioceptive signals influence this process is therefore critical for designing effective rehabilitation technologies. Systems such as Stimel-03 demonstrate how detecting weak neuromuscular activation, providing real-time feedback, and delivering synchronized stimulation can transform voluntary effort into meaningful therapeutic movement.

By integrating physiological monitoring, patient engagement, and targeted stimulation, intention-driven rehabilitation platforms support neuroplastic adaptation and functional recovery while establishing the sensory motor feedback loop that underpins modern neuromuscular rehabilitation.

---

## References

1. Biasucci A, Leeb R, Iturrate I, Perdakis S, Al-Khodairy A, Corbet T, et al. Brain-actuated functional electrical stimulation elicits lasting arm motor recovery after stroke. *Nat Commun.* 2018;9:2421.
2. Marquez-Chin C, Popovic MR. Functional electrical stimulation therapy for restoration of motor function after spinal cord injury and stroke: a review. *Biomed Eng Online.* 2020;19:34.

3. Eraifej J, Clark W, France B, Desando S, Moore D. Effectiveness of upper limb functional electrical stimulation after stroke: a systematic review and meta-analysis. *Clin Rehabil.* 2020;34(9):1188-1199.
4. Chuang LL, Chen YL, Chen CC, Li YC, Wong AM, Hsu AL, Chang YJ. Effect of EMG-triggered neuromuscular electrical stimulation with bilateral arm training on hemiplegic shoulder pain and arm function after stroke: a randomized controlled trial. *J Neuroeng Rehabil.* 2017 Nov 28;14(1):122.
5. Ambrosini E, Ferrante S, Pedrocchi A, Ferrigno G, Molteni F. Cycling induced by electrical stimulation improves motor recovery in postacute hemiparetic patients: a randomized controlled trial. *Neurorehabil Neural Repair.* 2020;34(6):505-516.
6. Knutson JS, Fu MJ, Sheffler LR, Chae J. Neuromuscular electrical stimulation for motor restoration in hemiplegia. *Phys Med Rehabil Clin N Am.* 2021;32(2):317-339.
7. Thrasher TA, Zivanovic V, McIlroy W, Popovic MR. Rehabilitation of reaching and grasping function in severe hemiplegic patients using functional electrical stimulation therapy. *Neurorehabil Neural Repair.* 2020;34(7):573-581.
8. Alon G, Levitt AF, McCarthy PA. Functional electrical stimulation enhancement of upper extremity functional recovery during stroke rehabilitation: a systematic review. *Am J Phys Med Rehabil.* 2021;100(11):1035-1044.