

STIMEL - 03 WHITE PAPER SERIES

---

White Paper 3

# Biofeedback, Proprioception, and Sensorimotor Reintegration

*The Role of Sensory Feedback in Motor Recovery*

---

## Abstract

---

Recovery of voluntary movement after neurological injury depends on restoring the sensorimotor feedback loops that connect the brain, muscles, and sensory receptors. After stroke and other neurological disorders, patients frequently retain motor intent that generates low amplitude EMG signals but insufficient force to produce visible movement. When these signals are detected and reinforced through biofeedback and stimulation assisted movement, the nervous system receives sensory confirmation that supports motor relearning and neuroplastic adaptation. This white paper explains the neurophysiology of proprioception, how neurological injury disrupts sensorimotor integration, and how biofeedback driven rehabilitation systems rebuild these pathways. It also describes how the Stimel-03 platform integrates EMG detection, adaptive signal processing, visual biofeedback, and electrical stimulation to convert weak voluntary activation into reinforced sensorimotor learning.


## Sensorimotor Control and the Role of Proprioception

---

Human movement depends on a continuous sensorimotor control loop linking the central nervous system and the musculoskeletal system. Motor commands generated in cortical and subcortical brain regions travel through descending neural pathways to activate specific muscles. When muscles contract, sensory receptors located in muscles, tendons, and joints transmit information back to the brain.

Three related sensory mechanisms contribute to this feedback process. Proprioception refers to the detection of limb position and muscle force through receptors such as muscle spindles and Golgi tendon organs. Kinesthesia refers specifically to the perception of limb motion. Together these signals allow the nervous system to evaluate whether the executed movement matches the intended motor command.

During motor learning the brain continuously compares predicted movement outcomes with incoming sensory feedback. When discrepancies occur, error based updating adjusts future motor commands. This predictive process gradually refines movement accuracy and coordination. Efficient motor control therefore depends on reliable sensory feedback within the sensorimotor loop [1,2].




**Figure 1. The Human Sensorimotor Control Loop** Motor command generated in the brain → neural transmission through descending pathways → muscle contraction → proprioceptive receptor activation → sensory feedback returning to the brain for comparison between intended and executed movement.

## Sensorimotor Disruption After Neurological Injury

Neurological injuries such as stroke disrupt both motor output and sensory integration. Damage to cortical or subcortical pathways weakens descending neural signals that normally activate muscles. Patients may still generate motor intent, but the resulting neuromuscular activation is often too weak to produce visible movement.

Without successful movement the nervous system receives limited proprioceptive feedback confirming that the intended action occurred. This absence of sensory confirmation prevents effective error based updating and slows reinforcement of motor pathways responsible for voluntary movement. Patients may repeatedly attempt movements without receiving the sensory signals necessary to strengthen the underlying neural circuits.

This disruption of the sensorimotor loop represents one of the major barriers to early functional recovery after stroke [5].



**Figure 2. Disrupted Sensorimotor Loop After Stroke** Motor intent → weak EMG activation → insufficient muscle contraction → reduced proprioceptive feedback → incomplete reinforcement of motor pathways.

## Biofeedback Driven Sensorimotor Reintegration

Biofeedback technologies restore elements of the sensorimotor loop by making physiological signals associated with voluntary activation visible and measurable. Surface electromyography detects electrical activity produced by motor unit activation, even when the contraction is too weak to generate visible movement.

Displaying these signals visually allows patients to recognize successful activation attempts and adjust their effort accordingly. This feedback transforms otherwise invisible neuromuscular signals into actionable information that guides motor learning.

When biofeedback is combined with electrical stimulation that assists the attempted movement, the system recreates the missing sensory confirmation required for effective motor learning. Repeated pairing of voluntary intent, assisted movement, and sensory feedback reinforces the neural pathways responsible for the action [6].

## Stimel-03 and Reinforcement of Motor Intent

Stimel-03 implements this neurophysiological framework through an integrated closed loop rehabilitation platform that links EMG detection, visual biofeedback, and stimulation assisted movement.

Surface electrodes detect EMG signals generated during voluntary activation attempts. These signals are processed through the BioRhythmIQ architecture which isolates physiologically meaningful motor unit activity while reducing environmental noise and motion artifacts.

The processed signal is evaluated against adaptive activation thresholds that adjust dynamically according to the patient's baseline activity and performance during therapy sessions. When EMG activity exceeds the activation threshold the system triggers electrical stimulation that assists the intended movement.

The assisted contraction activates proprioceptive receptors within muscles and tendons, generating sensory feedback that travels back to the central nervous system. This feedback confirms successful execution of the intended movement and reinforces the neural pathway responsible for the action.

Many conventional EMG triggered stimulation systems rely on fixed activation thresholds and relatively simple signal conditioning. These designs may fail to detect early voluntary activation signals when EMG amplitudes are extremely small, which is common during the early stages of neurological recovery. When thresholds are fixed, weak but meaningful activation attempts may not trigger stimulation, preventing reinforcement of emerging motor pathways.

Stimel-03 addresses this limitation through continuous signal evaluation and adaptive threshold logic. The system monitors the patient's EMG baseline during therapy and adjusts activation thresholds relative to the patient's current signal capability. This allows the platform to respond to genuine voluntary intent even when visible movement has not yet emerged, enabling earlier participation in active rehabilitation and reinforcing the link between effort and movement.

Synchronizing stimulation with voluntary effort preserves the temporal relationship between neural drive and sensory feedback. Maintaining this timing is critical because neuroplastic adaptation depends on consistent pairing between motor intent and sensory confirmation.

## **BioRhythmiQ Signal Processing Architecture**

---

Detecting early voluntary activation presents significant engineering challenges because EMG amplitudes during early rehabilitation may be extremely small and easily obscured by noise or motion artifacts.

The BioRhythmiQ architecture uses a structured signal processing pipeline to isolate meaningful neuromuscular activity. EMG signals are acquired from surface electrodes placed over the target muscle groups. The signals pass through band pass filtering stages designed to isolate the frequency range associated with motor unit action potentials while suppressing environmental interference.

Artifact reduction algorithms then remove signal distortions caused by electrode movement or external electrical sources. After filtering and artifact suppression the system evaluates EMG amplitude relative to patient specific baseline activity using adaptive threshold logic.

Threshold adaptation is based on continuous comparison between real time EMG activity and the patient's recent signal history. Rather than using a fixed trigger level, the system dynamically adjusts the activation threshold within a range that reflects the patient's current neuromuscular capability. This approach improves detection reliability across therapy sessions and reduces the likelihood that meaningful activation attempts will be missed.

Low latency signal processing preserves the temporal coupling between voluntary activation and stimulation delivery. When stimulation occurs during the patient's attempted contraction, proprioceptive feedback is

temporally aligned with the original motor command. This temporal alignment strengthens associative motor learning and improves reinforcement of the intended movement pathway.

This architecture enables detection of extremely weak voluntary activation signals that would otherwise remain clinically unnoticed during early rehabilitation.

This thresholding approach continuously adjusts to fluctuations in signal strength across therapy sessions. When EMG activity exceeds the adaptive threshold the trigger logic activates stimulation assistance. Low latency signal processing ensures that stimulation occurs during the patient's voluntary activation attempt, preserving the temporal pairing between neural intent and sensory feedback that underlies associative motor learning.

This architecture enables detection of extremely weak voluntary activation signals that would otherwise remain clinically unnoticed during early rehabilitation.

## Clinical Applications

---

Stimel-03 supports multiple rehabilitation scenarios in neurological therapy. In upper limb rehabilitation EMG guided stimulation assists wrist extension, finger opening, and grasp preparation during occupational therapy exercises. In lower limb rehabilitation the system assists ankle dorsiflexion and knee stabilization during gait training.

Because weak EMG signals frequently appear before visible movement returns, early detection allows patients to begin active rehabilitation sooner. Earlier activation increases the number of meaningful repetitions during therapy sessions and strengthens the sensorimotor learning process.

## Clinical Impact

---

Rehabilitation systems that combine biofeedback with stimulation assisted movement allow patients to actively participate in therapy rather than passively receiving stimulation. When voluntary activation attempts consistently produce assisted movement patients experience a direct connection between effort and outcome.

This connection increases engagement, improves therapy intensity, and reinforces the neural circuits responsible for voluntary movement. Repeated activation of these circuits promotes cortical reorganization and progressive functional recovery [8,9].

## Key Takeaways for Clinicians

---

- Weak voluntary EMG signals often appear before visible movement returns after neurological injury. Detecting these signals allows patients to begin active rehabilitation earlier.
- EMG driven biofeedback helps patients recognize successful activation attempts and improves engagement during therapy.

- Synchronizing electrical stimulation with voluntary intent restores the sensorimotor reinforcement cycle required for neuroplastic recovery.
- Stimel-03 detects low amplitude EMG signals and converts voluntary intent into assisted movement, enabling early stage rehabilitation training.
- Repeated pairing of voluntary activation, assisted movement, and sensory feedback strengthens neural pathways responsible for functional motor recovery.

## Conclusion

---

Motor recovery after neurological injury requires rebuilding the sensorimotor loops that connect the brain and the musculoskeletal system. Technologies that detect weak voluntary activation, provide visual feedback, and reinforce movement through stimulation restore the sensory confirmation required for effective motor learning.

By integrating EMG detection, adaptive signal processing, biofeedback visualization, and synchronized stimulation, Stimel-03 transforms the sensorimotor loop into a structured therapeutic process that supports neuroplastic adaptation and functional recovery.

## References

---

1. Proske U, Gandevia SC. The proprioceptive senses: their roles in signaling body shape, body position and movement. *Physiol Rev.* 2012.
2. Shadmehr R, Krakauer JW. A computational neuroanatomy for motor control. *Exp Brain Res.* 2008.
3. Goble DJ. Proprioceptive acuity assessment via joint position matching. *J Neurosci Methods.* 2010.
4. Aman JE, Elangovan N, Yeh IL, Konczak J. The effectiveness of proprioceptive training for improving motor function. *Front Hum Neurosci.* 2015.
5. Langhorne P, Bernhardt J, Kwakkel G. Stroke rehabilitation. *Lancet.* 2011.
6. Stanton R, Ada L, Dean CM, Preston E. Biofeedback improves performance in stroke rehabilitation. *J Physiother.* 2017.
7. Popovic MR, Kapadia N, Zivanovic V. Functional electrical stimulation therapy of voluntary grasping versus conventional rehabilitation. *Neurorehabil Neural Repair.* 2011.
8. Eraifej J, Clark W, France B. Effectiveness of upper limb functional electrical stimulation after stroke. *Clin Rehabil.* 2020.
9. Knutson JS, Fu MJ, Sheffler LR. Neuromuscular electrical stimulation for motor restoration in hemiplegia. *Phys Med Rehabil Clin N Am.* 2021.