

STIMEL - 03 WHITE PAPER SERIES

White Paper 6

Closed Loop Neurorehabilitation Systems

From Motor Intent to Functional Movement

Abstract

Effective neurorehabilitation depends on restoring the biological learning loop that links motor intention, movement execution, and sensory feedback. After neurological injury such as stroke, patients frequently attempt to move but cannot generate sufficient muscle contraction to produce visible movement or proprioceptive confirmation. When the brain does not receive sensory feedback confirming that an attempted action produced a result, reinforcement of the motor pathway weakens and recovery slows. Closed loop rehabilitation systems address this gap by detecting voluntary neural signals and converting them into assisted movement within a coherent time window. This paper explains the principles of closed loop neurorehabilitation and describes how the Stimel-03 platform operationalizes these mechanisms through EMG intent detection, BioRhythmIQ signal processing, adaptive activation thresholds, and synchronized stimulation.

Clinical Problem or Rehabilitation Gap

Motor recovery after neurological injury requires the brain to reconnect motor intention with physical action. In the early phases of stroke rehabilitation many patients generate genuine descending motor drive but cannot produce sufficient muscle activation to create visible movement. Surface electromyography often detects weak voluntary motor unit activity even when clinical observation suggests paralysis [3].

When attempted movements fail to produce sensory confirmation the brain receives little evidence that the motor command produced an effect. The absence of proprioceptive reinforcement disrupts the association between intention and movement that is required for neuroplastic adaptation. Rehabilitation systems must therefore restore the link between motor intent and sensory feedback so that attempted actions again produce meaningful reinforcement of motor circuits.

Scientific or Physiological Basis

Neuroplastic recovery depends on repeated activation of motor pathways within a coherent sensorimotor loop. When intention, movement, and sensory feedback occur together within a biologically meaningful time window, synaptic connections associated with the movement become stronger. This associative process underlies motor relearning following neurological injury.

Associative motor learning follows Hebbian reinforcement principles in which temporally aligned neural activity strengthens synaptic connections. When voluntary motor intent is immediately followed by movement and proprioceptive confirmation, the brain links the command with its sensory consequence. Preserving this temporal alignment is therefore critical for rehabilitation systems designed to promote functional recovery.

Surface electromyography provides a practical method for detecting voluntary motor unit activation. Voluntary EMG signals associated with motor unit action potentials typically occupy a frequency band of approximately 20–500 Hz, allowing physiologic activity to be distinguished from motion artifacts and environmental noise [3]. Even in early stroke recovery, small bursts of EMG activity can indicate genuine voluntary motor intent.

Detecting these weak signals is critical because they represent the earliest evidence that descending motor pathways remain partially functional.

Mechanism of Closed Loop Neurorehabilitation

Closed loop neurorehabilitation systems respond directly to the patient's voluntary neural activity. Instead of delivering stimulation according to a predetermined schedule, the system detects the patient's attempt to move and assists the intended action.




Figure 1. Conceptual Closed Loop Neurorehabilitation Cycle This cycle illustrates how closed loop systems reconnect the patient's intention with sensory confirmation, enabling associative motor learning even when voluntary contraction is initially weak.

The closed loop rehabilitation model involves four sequential stages.

- Detection of voluntary motor intent
- Signal interpretation through real time processing
- Conversion of intent into assisted contraction
- Sensory feedback confirming the movement

When these stages occur within a coherent time window the brain can associate intention with movement and strengthen the corresponding neural pathway. This alignment is essential for effective motor learning and rehabilitation.

Stimel-03 Implementation

Stimel-03 implements the closed loop rehabilitation model through an integrated system architecture combining EMG detection, BioRhythmIQ signal processing, adaptive activation thresholds, and synchronized stimulation.

Surface EMG acquisition captures voluntary motor unit activity generated by the patient's attempt to move.

BioRhythmIQ signal processing isolates physiologic EMG signals while suppressing interference from motion artifacts, electrode movement, muscle cross talk, and environmental electrical noise such as 50 or 60 Hz power line interference [2]. Through filtering and artifact suppression the system improves the signal to noise ratio and increases the reliability of voluntary intent detection.

Adaptive activation thresholds continuously adjust to the patient's current neuromuscular state [4]. Rather than relying on fixed trigger levels, the system evaluates recent EMG signal characteristics and updates the activation threshold dynamically. This allows the platform to respond to weak or fluctuating voluntary activation that conventional EMG triggered systems often miss.

When voluntary EMG activity crosses the adaptive threshold the system delivers stimulation that assists the intended movement. System latency remains within a short physiological window so that voluntary intent, assisted contraction, and proprioceptive feedback are perceived as part of the same motor event. Maintaining this temporal alignment strengthens neural reinforcement and supports effective motor relearning.

Stimel-03 Differentiation from Conventional EMG Triggered FES

Conventional EMG triggered stimulation systems typically treat EMG as a binary trigger event based on a fixed threshold level. Weak voluntary signals that fall below this threshold are often ignored, particularly during early stage recovery when motor unit recruitment is minimal.

Stimel-03 interprets EMG as a dynamic signal representing evolving motor intent. Adaptive thresholds adjust continuously to the patient's capability, allowing detection of low amplitude voluntary activity. BioRhythmIQ signal processing improves discrimination between true motor unit activation and noise, while low latency stimulation preserves the biological timing required for associative learning. Together these capabilities allow intention driven rehabilitation even when visible movement is absent.

Clinical Applications

Several rehabilitation scenarios illustrate the clinical use of closed loop neurorehabilitation.

Early stroke wrist extension

A patient attempts wrist extension but produces only weak EMG activity without visible movement. Stimel-03 detects the motor intent and triggers stimulation, allowing the patient to see and feel the intended movement. Repeated intention driven cycles reinforce voluntary control of the extensor muscles.

Hand opening during upper limb retraining

During grasp release training a patient generates low amplitude voluntary activation that is insufficient to open the hand. The system assists the movement and reinforces the connection between voluntary effort and sensory confirmation.

Ankle dorsiflexion during gait rehabilitation

A patient attempting dorsiflexion during walking produces minimal activation of the tibialis anterior muscle. Stimel-03 converts the detected effort into assisted contraction, improving foot clearance and reinforcing the motor pattern required for gait.

Clinical Scenario Example

A 64 year old patient 10 days after ischemic stroke attempts wrist extension during therapy but produces no visible movement. Surface EMG reveals microvolitional bursts of approximately 12–20 μV in the extensor muscle group. Stimel-03 detects this early voluntary activity and delivers synchronized stimulation that produces visible wrist extension. The patient observes the movement and experiences proprioceptive feedback confirming the intended action. Repeated intention driven cycles strengthen the association between voluntary effort and motor output and support the gradual return of independent wrist extension.

Clinical Impact

Closed loop rehabilitation systems allow therapy to begin earlier in recovery when voluntary EMG activity exists but visible movement has not yet returned [1]. By converting weak voluntary intent into assisted movement these systems restore the biological learning loop that supports neuroplastic adaptation.

Patients experience a stronger connection between attempted actions and observed outcomes, which increases engagement and encourages repeated activation of recovering motor circuits. Increased repetition and consistent sensorimotor reinforcement are key drivers of functional recovery in neurorehabilitation.

Integration into therapist workflow is straightforward. Surface electrodes are applied in the same manner as standard EMG or stimulation therapy, and sessions can be incorporated into conventional physiotherapy or occupational therapy exercises. Therapists guide task oriented movements while the system detects voluntary intent and provides synchronized assistance, allowing closed loop reinforcement to occur during normal rehabilitation activities.

Key Takeaways for Clinicians

- Motor recovery requires restoration of the intention movement feedback learning loop.
- Weak voluntary EMG signals may exist even when visible movement is absent.
- Closed loop rehabilitation systems convert detected intent into assisted movement.
- Stimel-03 operationalizes closed loop neurorehabilitation through EMG detection, BioRhythmIQ signal processing, adaptive thresholds, and synchronized stimulation.
- These mechanisms allow patients to participate actively in rehabilitation earlier in recovery.

Conclusion

Closed loop neurorehabilitation restores the fundamental biological learning cycle required for recovery of voluntary motor control after neurological injury. When patients attempt to move but cannot generate visible contraction, the brain receives little sensory confirmation that the motor command produced a result. Without this reinforcement, opportunities for neuroplastic adaptation are limited and motor recovery may slow.

Closed loop rehabilitation systems address this problem by detecting voluntary motor intent and converting it into assisted movement that occurs within a coherent physiological time window. By synchronizing stimulation with the patient's own neuromuscular activation, these systems reconnect intention, movement execution, and sensory feedback. This temporal alignment allows the nervous system to associate voluntary effort with successful movement and reinforces the neural circuits responsible for motor control.

Stimel-03 operationalizes this rehabilitation model through integrated EMG detection, BioRhythmIQ signal processing, adaptive activation thresholds, and synchronized stimulation delivery. These mechanisms enable the system to detect weak voluntary neuromuscular signals that often appear early in recovery and translate them into meaningful therapeutic movement. By restoring the sensorimotor learning loop and enabling earlier participation in active rehabilitation, closed loop systems such as Stimel-03 support progressive strengthening of voluntary motor pathways and contribute to functional recovery following neurological injury.

References

1. Khan MA, Fares H, Ghayvat H, et al. A systematic review on functional electrical stimulation based rehabilitation systems for upper limb post stroke recovery. *Front Neurol.* 2023;14:1272992.
2. Jung H, Lee J, Kim H. Volitional EMG estimation method during functional electrical stimulation by dual-channel surface EMGs. *Sensors.* 2021;21(23):8015.
3. Lv W, Liu K, Zhou P, Huang F and Lu Z (2023) Surface EMG analysis of weakness distribution in upper limb muscles post-stroke. *Front. Neurol.* 14:1135564. doi: 10.3389/fneur.2023.1135564
4. Carvalho CR, Fernández JM, Del-Ama AJ, Barroso FO, Moreno JC. Review of electromyography onset detection methods for real-time control of robotic exoskeletons. *Journal of NeuroEngineering and Rehabilitation.* 2023.