

STIMEL - 03 WHITE PAPER SERIES

White Paper 7

Detecting Motor Intent in Early Neurorehabilitation

Surface EMG as a Gateway to Intention-Driven Therapy

Abstract

Early neurological recovery frequently begins with weak voluntary neural signals that are insufficient to generate visible movement. These microvolitional signals represent genuine motor intent and indicate that descending motor pathways may remain partially functional even when clinical observation suggests paralysis. Conventional rehabilitation approaches often fail to detect these signals, delaying the initiation of intention-driven therapy. Surface electromyography provides a sensitive method for identifying early voluntary motor unit activity during attempted movement. This paper explains the neurophysiological basis of weak voluntary EMG signals, the challenges associated with detecting them, and how rehabilitation systems such as Stimel-03 convert detected motor intent into assisted movement that restores the sensorimotor learning loop.

Introduction: Clinical Problem in Early Motor Recovery

Patients recovering from stroke or other neurological injury frequently attempt to move before they can produce visible muscle contraction. In these early stages voluntary neural signals may still reach spinal motor neurons, but the resulting muscle activation is too weak to generate functional movement.

Clinically this often appears as complete paralysis. Neurophysiological studies, however, demonstrate that voluntary motor unit activity can still be detected in many patients using surface electromyography [1]. These signals indicate that the brain is attempting to activate the affected muscle even when contraction is not visible.

When weak voluntary signals go undetected, patients receive no sensory confirmation that their movement attempt produced an effect. The nervous system therefore receives fewer opportunities to reinforce the intended motor action. Detecting early voluntary motor signals is therefore essential for initiating rehabilitation during the most plastic phase of recovery.

Scientific Basis of Weak Voluntary EMG Signals

Motor recovery after neurological injury rarely occurs as a sudden return of full movement. Instead, recovery typically begins with partial activation of surviving neural pathways.

After stroke, partial corticospinal tract conduction may still transmit descending motor commands from the cortex to spinal motor neurons. However, the number of activated motor units is often insufficient to generate visible contraction. In this situation central motor drive is present but peripheral expression of movement remains weak. Failure of adequate motor-unit summation means that muscle activation may only appear as small electrical discharges detectable with EMG rather than visible motion.

During early recovery several physiological conditions may therefore be present:

- Descending cortical signals still reach spinal motor neurons
- Motor unit recruitment is sparse and inconsistent
- Muscle contraction remains too weak to produce visible movement

Surface electromyography can detect the electrical activity produced by these early motor-unit activations even when contraction remains below the threshold required for visible movement. These low-amplitude signals represent genuine voluntary effort and therefore provide an objective indicator of motor intent.

Detecting these signals allows rehabilitation systems to respond directly to the patient's attempted movement rather than relying on externally triggered stimulation.

Recent clinical studies have also shown that surface EMG can reveal neuromotor activity patterns during early stroke recovery that are not detectable during routine clinical examination, reinforcing the value of EMG as an early biomarker of voluntary motor activation [4].

Characteristics of Weak Voluntary EMG Signals

- Weak voluntary EMG signals recorded during early neurorehabilitation typically show several identifiable characteristics.
- Low amplitude. Signal amplitudes are often extremely small, commonly in the range of approximately 10–30 microvolts.
- Irregular activation patterns. Motor unit firing may be intermittent and unstable during early recovery.
- High susceptibility to noise. Because the signals are small they can be easily obscured by motion artifacts, electrode impedance fluctuations, or environmental electrical interference.
- Reliable detection therefore requires signal processing approaches capable of distinguishing weak voluntary activity from background noise.

Detecting Motor Intent Using Surface EMG

Surface EMG provides a practical method for detecting voluntary muscle activation in clinical rehabilitation environments. Voluntary EMG signals produced by motor unit activity typically occupy a frequency band between approximately 20 and 500 Hz [2]. Within this band signal features such as amplitude, timing, and activation patterns can be analyzed to identify voluntary intent.

Recent EMG-intention detection studies further demonstrate that advanced signal-processing methods can reliably identify weak voluntary activation patterns associated with attempted movement, enabling intention-driven rehabilitation systems to respond in real time to patient effort [5].

When rehabilitation systems detect these signals during attempted movement, therapy can begin even before visible movement returns. Detecting early motor intent enables the conversion of voluntary neural activity into assisted movement, restoring the connection between intention and sensory feedback.

Stimel-03 Detection and Reinforcement Mechanism

Stimel-03 integrates surface EMG detection with BioRhythmIQ, the system's adaptive EMG interpretation engine, to identify weak voluntary motor unit activity during early rehabilitation.

Surface electrodes capture neuromuscular electrical activity generated when the patient attempts to activate the affected muscle. BioRhythmIQ signal processing analyzes the signal continuously, isolating physiologic EMG activity while suppressing motion artifacts, cross-talk between adjacent muscles, and environmental electrical interference. The system evaluates signal trends and baseline variability in real time, enabling interpretation of weak voluntary activity patterns rather than relying solely on simple amplitude thresholds.

Adaptive activation thresholds continuously adjust to the patient's real-time neuromuscular state. Unlike conventional EMG-triggered stimulation systems that rely on fixed trigger thresholds, Stimel-03 interprets EMG as a dynamic physiological signal. This allows the system to detect micro-volitional motor unit activity that static threshold systems frequently miss.

When voluntary EMG activity crosses the adaptive threshold, stimulation is delivered to assist the intended movement. The resulting contraction produces visible movement and generates proprioceptive feedback confirming that the attempted action produced an effect.

Equally important is temporal alignment. Once voluntary EMG activity is detected, stimulation delivery occurs within a biologically meaningful reinforcement window typically well below 150 milliseconds. Preserving this timing relationship allows the brain to associate voluntary effort with the resulting movement and sensory feedback, strengthening the neural pathways involved in motor control through associative neuroplastic mechanisms [3].

Differentiation from Conventional EMG-Triggered Stimulation

Many conventional EMG-triggered stimulation systems treat EMG as a simple binary threshold event. When the signal exceeds a preset amplitude stimulation is delivered, and when it does not the attempt is ignored. In early rehabilitation weak voluntary signals often remain below these fixed trigger thresholds, meaning genuine motor attempts may fail to initiate therapy.

Stimel-03 approaches EMG detection differently. BioRhythmIQ signal processing evaluates the EMG signal continuously and adapts activation thresholds based on the patient's moment-to-moment signal characteristics. This dynamic interpretation enables detection of weak voluntary activity that conventional systems frequently overlook. By converting these early intention signals into synchronized assisted movement, Stimel-03 restores the functional relationship between motor intent, movement execution, and sensory confirmation.

Clinical Scenario Example

A patient two weeks after stroke attempts wrist extension but produces no visible movement. Surface EMG reveals weak voluntary bursts of approximately 15 microvolts in the wrist extensor muscles. Stimel-03 detects these signals and triggers stimulation that produces visible wrist extension. The patient observes the movement and experiences proprioceptive feedback confirming the intended action. Repeated cycles of

intention-driven stimulation reinforce the association between voluntary effort and movement and support progressive recovery of independent motor control.

Clinical Impact

- Detecting weak voluntary EMG signals has several important implications for neurorehabilitation.
- Therapy can begin earlier in recovery when voluntary neural activity first appears.
- Patients receive immediate confirmation that their movement attempt produces a visible effect.
- Each successful intention-driven movement reinforces the neural pathways associated with the intended action.

These mechanisms increase meaningful repetitions of motor activity, which is a central driver of neuroplastic recovery.

Key Takeaways for Clinicians

- Weak voluntary EMG signals may be present even when visible movement is absent.
- Surface EMG provides a sensitive method for detecting early motor intent.
- Detecting voluntary EMG allows rehabilitation systems to respond directly to the patient's attempted movement.
- Stimel-03 converts detected motor intent into assisted movement, restoring the intention-movement-feedback learning loop.
- Intention-driven therapy enables rehabilitation to begin earlier during the most plastic phase of neurological recovery.

Conclusion

Early detection of voluntary neural activity reshapes how rehabilitation can begin. By identifying weak EMG signals that represent genuine motor intent, clinicians can initiate intention-driven therapy even before visible movement returns. Systems such as Stimel-03 translate these early neural signals into assisted movement, reconnecting motor intention with sensory feedback. This restoration of the intention-movement-feedback loop enables rehabilitation to begin earlier, increases meaningful repetitions, and strengthens the neural pathways required for functional recovery.

References

1. Steele KM, Lee SSM, Grudzińska K, Apkarian A, Brinegar K, Wolf SL, et al. Muscle activity after stroke: perspectives on deploying surface electromyography in acute care. *Front Neurol.* 2020;11:576757.
2. Jung H, Lee J, Kim H. Volitional EMG estimation method during functional electrical stimulation by dual-channel surface EMGs. *Sensors.* 2021;21(23):8015.
3. Kleim JA, Jones TA. Principles of experience-dependent neural plasticity: implications for rehabilitation after brain damage. *J Speech Lang Hear Res.* 2008;51:S225-S239.

4. Lv W, Liu K, Zhou P, Huang F and Lu Z (2023) Surface EMG analysis of weakness distribution in upper limb muscles post-stroke. *Front. Neurol.* 14:1135564. doi: 10.3389/fneur.2023.1135564
5. Carvalho CR, Fernández JM, Del-Ama AJ, Barroso FO, Moreno JC. Review of electromyography onset detection methods for real-time control of robotic exoskeletons. *Journal of NeuroEngineering and Rehabilitation.* 2023.